

Product datasheet for **AP55142SU-N**

Peptidyl-prolyl cis-trans isomerase A (PPIA/PPIase 1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/500-1/5000. Immunohistochemistry: 1/100-1/500.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from an internal domain of Human PPIase1
Specificity:	This antibody reacts with Human 18 kDa protein. Cross reacts with proteins from many other species due to the homology sequences.
Formulation:	State: Serum State: Lyophilized serum Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store the antibody at -20°C. Store the reconstituted antibody at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<u>Q567Q0</u>



[View online »](#)