

## Product datasheet for **AP55151SU-N**

### Phosphoglycerate kinase (PGK) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot:</b> 1/200-1/2000. <b>Immunohistochemistry:</b> 1/50-1/500.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of PGK from <i>Trypanasoma brucei</i>
Specificity:	Reacts with <i>Trypanasoma Brucei</i> 101 kDa protein.
Formulation:	State: Serum State: Lyophilized powder Preservative: None
Reconstitution Method:	Restore in distilled water
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store the antibody at -20°C. Store reconstituted antibody at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<u><a href="#">P07377</a></u>



[View online »](#)