

Product datasheet for **AP55153SU-N**

plc-1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot: 1/500-1/2000.
Reactivity:	C. elegans
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from an internal domain of Plc-1 protein
Specificity:	Reacts with 212 kDa Plc-1 protein from <i>C. elegans</i> .
Formulation:	State: Serum State: Lyophilized powder Preservative: None
Conjugation:	Unconjugated
Storage:	Store lyophilized antibody at -20°C Store reconstituted antibody at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	G5EF18



[View online »](#)