

Product datasheet for **AP55154SU-N**

plcd Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Sea Urchin
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal parts of PLC delta
Specificity:	Reacts with purple sea urchin PLC delta 85kDa protein.
Formulation:	State: Serum State: Lyophilized powder Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Lyophilized power stable for a minimum of 2 years at -20°C. Store reconstituted antibody at 2-8°C or (in aliquots) at -20°C for longer.
Stability:	Shelf life: Six months from despatch.
Database Link:	Q5QIB9



[View online »](#)