

## Product datasheet for **AP55155SU-N**

### plc-3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot:</b> 1/500 - 1/2000.
Reactivity:	C. elegans
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of the Plc-3 protein
Specificity:	This antibody reacts with 154 kDa Plc-3 protein from C. elegans.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore with distilled water.
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)