

Product datasheet for **AP55156SU-N**

plc-4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	Western Blot: 1/500-1/2000.
Reactivity:	<i>C. elegans</i>
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from an internal domain of Plc-4 protein
Specificity:	This antibody reacts with 85 kDa Plc-4 protein from <i>C. elegans</i> .
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	Q21754



[View online »](#)