

Product datasheet for **AP55165SU-N**

PABPC1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections: 1/100 - 1/1000. Western Blot: 1/500 - 1/5000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of human PABP-1 protein
Specificity:	This antibody reacts with 70 kDa PABP-1 protein.
Formulation:	State: Serum State: Lyophilized Serum
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	poly(A) binding protein, cytoplasmic 1
Database Link:	Entrez Gene 18458 Mouse P29341
Synonyms:	PAB1, PABP-1, PABP1, Polyadenylate-binding protein 1, PABPC2



[View online »](#)