

Product datasheet for **AP55166SU-N**

PABPC4 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections: 1/100 - 1/1000. Western Blot: 1/500 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of human PABP-4 protein
Specificity:	This antibody reacts with 70 kDaPABP-4 protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	poly(A) binding protein cytoplasmic 4
Database Link:	Entrez Gene 8761 Human Q13310
Synonyms:	PABP4, PABP-4, Polyadenylate-binding protein 4, APP1, APP-1, iPABP



[View online »](#)