

## Product datasheet for **AP55183PU-N**

### KCNE1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry:</b> 1/100 - 1/1000. <b>Western Blot:</b> 1/500 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of KCNE1 protein
Specificity:	This antibody reacts with 14 kDa KCNE1 protein.
Formulation:	0.1M Tris, 0.1M Glycine, 2% Sucrose State: Purified State: Lyophilized purified powder
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Purification:	Affinity chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	potassium voltage-gated channel subfamily E regulatory subunit 1
Database Link:	<a href="#">Entrez Gene 3753 Human P15382</a>
Synonyms:	Minimal potassium channel, Delayed rectifier potassium channel subunit Isk



[View online »](#)