

## Product datasheet for **AP55194SU-N**

### PRNP Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry:</b> 1/200 - 1/1000. <b>Western Blot:</b> 1/200 - 1/1000.
Reactivity:	Bovine, Human, Rat, Sheep
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A 15 residue synthetic peptide derived from C-terminal domain of bovin PrP protein
Specificity:	This antibody reacts with full sequence PrP protein of brain tissue.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<a href="#">P23907</a>
Synonyms:	Major prion protein, PrP27-30, PrP33-35C, ASCR, PRNP, PRIP



[View online »](#)