

## Product datasheet for **AP55195PU-N**

### PRNP Rabbit Polyclonal Antibody

#### Product data:

|                       |  |
|-----------------------|--|
| Product Type:         | Primary Antibodies   |
| Applications:         | ELISA, IHC, WB   |
| Recommended Dilution: | <b>ELISA.</b><br><b>Immunohistochemistry:</b> 1/20 - 1/100.<br><b>Western Blot:</b> 1/20 - 1/100.  |
| Reactivity:           | Bovine, Human, Rat, Sheep  |
| Host:                 | Rabbit   |
| Clonality:            | Polyclonal   |
| Immunogen:            | A 15 residue synthetic peptide derived from C-terminal domain of bovin PrP protein   |
| Specificity:          | This antibody reacts with full sequence PrP protein of brain tissue.   |
| Formulation:          | Tris 0.1 M, glycine 0.1 M, Glycerol 50%, (pH 7.8)<br>State: Aff - Purified<br>State: Liquid purified Ig fraction<br>Preservative: 0.02% sodium azide |
| Purification:         | Affinity chromatography  |
| Conjugation:          | Unconjugated   |
| Storage:              | Store the antibody at -20°C.<br>Avoid repeated freezing and thawing.   |
| Stability:            | Shelf life: one year from despatch.  |
| Database Link:        | <a href="#">P23907</a>   |
| Synonyms:             | Major prion protein, PrP27-30, PrP33-35C, ASCR, PRNP, PRIP   |



[View online »](#)