

Product datasheet for **AP55199SU-N**

Senataxin (SETX) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot. Immunohistochemistry.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptides derived from internal domain of the human SETX protein
Specificity:	This antibody reacts with 303 kDa SETX protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	senataxin
Database Link:	Entrez Gene 23064 Human Q7Z333
Synonyms:	SEN1 homolog, SCAR1, ALS4, KIAA0625



[View online »](#)