

## Product datasheet for **AP55207SU-N**

### PSMA4 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry:</b> 1/50 - 1/500. <b>Western Blot:</b> 1/500 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of human PSA4
Specificity:	This antibody reacts with 29 kDa PSA4 protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	proteasome subunit alpha 4
Database Link:	<a href="#">Entrez Gene 5685 Human P25789</a>
Synonyms:	Proteasome subunit alpha type-4, Macropain subunit C9, Proteasome component C9, HC9, PSC9, Proteasome subunit L



[View online »](#)