

Product datasheet for **AP55213PU-N**

LYRIC (MTDH) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry: 1/50 - 1/200. Western Blot: 1/100 to 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of human LYRIC protein
Specificity:	This antibody reacts with 63 kDa protein LYRIC.
Formulation:	0.1M Tris, 0.1M Glycine, 2% Sucrose State: Purified State: Lyophilized purified powder
Reconstitution Method:	Restore in distilled water.
Concentration:	lot specific
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	metadherin
Database Link:	Entrez Gene 92140 Human Q86UE4
Synonyms:	3D3/LYRIC, Metastasis adhesion protein, AEG1, LYRIC, Metadherin, AEG-1



[View online »](#)