

Product datasheet for **AP55215SU-N**

NDRG3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry: 1/200 - 1/1000. Western Blot: 1/500 to 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of the human NDRG3 protein.
Specificity:	This antibody reacts with 41 kDa NDRG3 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	NDRG family member 3
Database Link:	Entrez Gene 57446 Human Q9UGV2
Synonyms:	FLJ13556



[View online »](#)