

Product datasheet for **AP55224SU-N**

Ptch1 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry: 1/50 - 1/500. Western Blot: 1/500 - 1/2000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of mouse PTC1
Specificity:	This antibody reacts with 160 kDa PTC1 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	patched 1
Database Link:	Entrez Gene 19206 Mouse Q61115
Synonyms:	PTCH, PTC1, PTC, Patched 1



[View online »](#)