

## Product datasheet for **AP55227SU-N**

### PPP1R15B Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Immunohistochemistry:</b> 1/50 - 1/200. <b>Western Blot:</b> 1/100 - 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of human PR15B protein
Specificity:	This antibody reacts with 79 kDa PR15B protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	protein phosphatase 1 regulatory subunit 15B
Database Link:	<u><a href="#">Entrez Gene 84919 Human</a></u> <u><a href="#">Q5SWA1</a></u>
Synonyms:	FLJ14744



[View online »](#)