

Product datasheet for **AP55228SU-N**

PTC2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/200-1/2000. Immunohistochemistry: 1/50-1/500.
Reactivity:	Yeast
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from the C-terminal domain of PP2C2
Specificity:	Reacts with Yeast 50 kDa Ptc2 protein
Formulation:	State: Serum State: Lyophilized serum Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	P39966



[View online »](#)