

## Product datasheet for **AP55232SU-N**

### SEC31B Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>Immunohistochemistry:</b> 1/100 - 1/500. <b>ELISA.</b> <b>Western Blot:</b> 1/500 - 1/5000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of human Sec31B
Specificity:	This antibody reacts with 128 kDa Sec31B protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	SEC31 homolog B, COPII coat complex component
Database Link:	<a href="#">Entrez Gene 25956 Human Q9NQW1</a>
Synonyms:	SEC31-related protein B, SEC31B-1, SEC31-like protein 2, SEC31L2



[View online »](#)