

## Product datasheet for **AP55235SU-N**

### Lubricin (PRG4) Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>Immunohistochemistry:</b> 1/50 - 1/200. <b>ELISA.</b> <b>Western Blot:</b> 1/200 - 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from C-term domain of PRG4 protein
Specificity:	This antibody reacts with 151 kDa PRG4 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	proteoglycan 4
Database Link:	<a href="#">Entrez Gene 10216 Human Q92954</a>
Synonyms:	Lubricin, PRG4, MSF, SZP



[View online »](#)