

Product datasheet for **AP55237SU-N**

FES Rabbit Polyclonal Antibody

Product data:

| | |
|------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, IHC, WB |
| Recommended Dilution: | ELISA. Immunohistochemistry: 1/100 - 1/1000. Western Blot: 1/1000 - 1/5000. |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic peptide derived from human Fes protein |
| Specificity: | This antibody reacts with 93 kDa Fes protein. |
| Formulation: | State: Serum State: Lyophilized serum |
| Reconstitution Method: | Restore in distilled water. |
| Conjugation: | Unconjugated |
| Storage: | Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |
| Gene Name: | FES proto-oncogene, tyrosine kinase |
| Database Link: | Entrez Gene 2242 Human P07332 |
| Synonyms: | FES, FPS, Proto-oncogene c-Fps |



[View online »](#)