

Product datasheet for **AP55239SU-N**

FLJ14213 (PRR5L) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/200. ELISA. Western Blot: 1/100 - 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of human Ptor-2 protein
Specificity:	This antibody reacts with 40 kDa Ptor-2 protein.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	proline rich 5 like
Database Link:	Entrez Gene 79899 Human Q6MZQ0
Synonyms:	PROTOR2, Proline-rich protein 5-like, FLJ14213



[View online »](#)