

Product datasheet for **AP55246SU-N**

RTN3 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. ELISA. Western Blot: 1/500 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of human RTN-3
Specificity:	This antibody reacts with 112 kDa RTN-3 protein. Cross reacts with human 2,3,4 & 5 isoforms.
Formulation:	0.1M Tris, 0.1M Glycine, 2% Sucrose State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	reticulon 3
Database Link:	Entrez Gene 10313 Human O95197
Synonyms:	ASYIP, NSPL2, NSPLII, NSP-like protein 2



[View online »](#)