

Product datasheet for **AP55248SU-N**

Retinol dehydrogenase 16 (RDH16) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. ELISA. Western Blot: 1/200 to 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from an internal domain of human RDH16
Specificity:	This antibody reacts with 36 kDa RDH16 protein.
Formulation:	0.1M Tris, 0.1M Glycine and 2% Sucrose State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	retinol dehydrogenase 16 (all-trans)
Database Link:	Entrez Gene 8608 Human O75452
Synonyms:	Retinol dehydrogenase 16, RoDH-4



[View online »](#)