

Product datasheet for **AP55262SU-N**

Securin (PTTG1) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Western Blot: 1/200 - 1/1000. Immunohistochemistry: 1/50 - 1/200.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of the human Securin protein
Specificity:	This antibody reacts with 22 kDa Securin protein.
Formulation:	0.1M Tris, 0.1M Glycine and 2% Sucrose State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Lyophilized power stable for a minimum of 2 years at -20°C. Store reconstituted antibodies at +4°C. For extended periods store in aliquots at -20°C. Antibodies are guaranteed for 6 month from date of receipt.
Gene Name:	pituitary tumor-transforming 1
Database Link:	Entrez Gene 9232 Human O95997
Synonyms:	PTTG, EAP1, TUTR1



[View online »](#)