

## Product datasheet for **AP55291SU-N**

### Slc12a7 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IF, IHC, WB
Recommended Dilution:	<b>Western Blot.</b> <b>Immunohistochemistry.</b> <b>Immunofluorescence.</b> <b>ELISA.</b>
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from internal domain of mouse Slc12A7 protein
Specificity:	This antibody reacts with Slc12a7 112kDa protein.
Formulation:	0.1M Tris, 0.1M Glycine and 2% Sucrose State: Serum State: Lyophilized powder
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	solute carrier family 12, member 7
Database Link:	<a href="#">Entrez Gene 20499 Mouse</a> <a href="#">Q9WVL3</a>
Synonyms:	KCC4



[View online »](#)