

## Product datasheet for **AP55294SU-N**

### Spermine Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	<b>ELISA.</b> <b>Western Blot:</b> 1/1200-1/5000. <b>Immunohistochemistry:</b> 1/50-1/500.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	KLH-Spermine
Specificity:	This antibody reacts with Spermine and cross reacts with Spermidine.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore with 200 µl distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)