

## Product datasheet for **AP55300SU-N**

### Spns3 pSer3 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot: 1/500 - 1/2000.
Reactivity:	Mouse
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from mouse C terminal region of Spinster-3 protein
Specificity:	This antibody reacts with 55 kDa SPNS3 protein.
Formulation:	0.1M Tris, 0.1M Glycine and 2% Sucrose State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	spinster homolog 3
Database Link:	<a href="#">Entrez Gene 77577 Mouse Q9D232</a>
Synonyms:	Protein spinster homolog 3



[View online »](#)