

Product datasheet for **AP55311SU-N**

Syntaxin 3 (STX3) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry. Immunofluorescence: 1/250. ELISA. Western Blot: 1/1000 - 1/10000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Syntaxin-3-GST fusion protein
Specificity:	This antibody reacts with 35 kDa STX3 protein.
Formulation:	0.1M Tris, 0.1M Glycine and 2% Sucrose State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	syntaxin 3
Database Link:	Entrez Gene 6809 Human Q13277
Synonyms:	Syntaxin-3, STX3A



[View online »](#)