

Product datasheet for **AP55313SU-N**

MAL2 Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/200. ELISA. Western Blot: 1/200 - 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from N-term domain of MAL2 protein
Specificity:	Reacts with human 19 kDa MAL2 protein. Cross reacts with mouse, rat and other species.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	mal, T-cell differentiation protein 2 (gene/pseudogene)
Database Link:	Entrez Gene 114569 Human Q969L2
Synonyms:	Protein MAL2, MAL-2



[View online »](#)