

Product datasheet for **AP55338SU-N**

Transglutaminase 4 (TGM4) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/200 - 1/500. ELISA. Western Blot: 1/200 -o 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	A 15 residue synthetic peptide derived from specific domain of Tgase 4 protein
Specificity:	This antibody reacts with TGase-4.
Formulation:	State: Serum State: Lyophilized serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	transglutaminase 4
Database Link:	Entrez Gene 7047 Human P49221
Synonyms:	TGase-4, TGP, Fibrinolygase, Prostate transglutaminase



[View online »](#)