

Product datasheet for **AP55369PU-N**

Pk (V5) Epitope Tag (GKPIP NPLLGLDST) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, WB
Recommended Dilution:	ELISA. Western Blot.
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	V5 synthetic peptide conjugated to KLH
Specificity:	This antibody reacts with V5-tagged fusion proteins.
Formulation:	0.1M Tris, 0.1M Glycine and 2% Sucrose State: Purified State: Lyophilized purified powder Preservative: None
Reconstitution Method:	Restore in distilled water.
Purification:	Affinity Chromatography on Protein A
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store the antibody at -20°C. Store reconstituted antibody at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Background:	The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 tag is usually used with all 14 amino acids (GKPIP NPLLGLDST), although it has also been used with a shorter 9 amino acid sequence (IPNPLLGLD).



[View online »](#)