

## Product datasheet for **AP55396SU-N**

### WDHD1 Rabbit Polyclonal Antibody

#### Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	Immunohistochemistry: 1/50 - 1/500. Western Blot: 1/500 - 1/2000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Nter domain of human WDHD1
Specificity:	This antibody reacts with human 125 kDa WDHD1 protein. Cross reacts with mouse protein.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	WD repeat and HMG-box DNA binding protein 1
Database Link:	<a href="#">Entrez Gene 11169 Human</a> <a href="#">O75717</a>
Synonyms:	AND1, AND-1



[View online »](#)