

Product datasheet for **AP55399SU-N**

DFNB31 (WHRN) Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry on frozen sections: 1/50 - 1/200. Western Blot: 1/100 - 1/1000.
Reactivity:	Human
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Synthetic peptide derived from Cter domain of human WHRN protein
Specificity:	This antibody reacts with 96 kDa WHRN protein. Cross reacts with mouse, rat and all protein isoforms.
Formulation:	State: Serum State: Lyophilized Serum
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Gene Name:	whirlin
Database Link:	Entrez Gene 25861 Human Q9P202
Synonyms:	WHRN, DFNB31, KIAA1526



[View online »](#)