

Product datasheet for **AP55410SU-N**

v-Jun Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	ELISA, EMSA, IHC, WB
Recommended Dilution:	ELISA. Immunohistochemistry on Frozen Sections: 1/50-1/500. Western Blot: 1/200-1/5000. Electrophoretic mobility shift assay (EMSA).
Reactivity:	Virus
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	v-Jun-GST fusion protein
Specificity:	This antibody reacts with avian and mammalian E.F (embryonic fibroblast) cells 32 kDa v-jun protein. Cross reacts with avian c-jun protein.
Formulation:	State: Serum State: Lyophilized Serum Preservative: None
Reconstitution Method:	Restore in distilled water.
Conjugation:	Unconjugated
Storage:	Store lyophilized at 2-8°C for 6 months or at -20°C long term. After reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C long term. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Database Link:	<u>P05411</u>



[View online »](#)