

Product datasheet for **AP60016PU-N**

PSAP Rabbit Polyclonal Antibody

Product data:

| | |
|------------------------|---|
| Product Type: | Primary Antibodies |
| Applications: | ELISA, IHC, WB |
| Recommended Dilution: | Immunohistochemistry: 1/200 - 1/1000. ELISA. Western Blot: 1/500 - 1/5000. |
| Reactivity: | Human |
| Host: | Rabbit |
| Clonality: | Polyclonal |
| Immunogen: | Synthetic peptide derived from an internal domain (Saposin-B) of human PSAP |
| Specificity: | This antibody reacts with 58 kDa PSAP protein. |
| Formulation: | 0.1M Tris, 0.1M Glycine and 2% Sucrose State: Purified State: Lyophilized purified powder |
| Reconstitution Method: | Restore in distilled water. |
| Concentration: | lot specific |
| Purification: | Affinity chromatography on Protein A |
| Conjugation: | Unconjugated |
| Storage: | Prior to reconstitution store at -20°C. Following reconstitution store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |
| Gene Name: | prosaposin |
| Database Link: | Entrez Gene 5660 Human P07602 |
| Synonyms: | Prosaposin, SAP1, GLBA |



[View online »](#)