

Product datasheet for **AR00019SU-N**

Borrelia afzelii Human Protein

Product data:

Product Type:	Native Proteins
Description:	Borrelia afzelii protein, 1 mg
Species:	Human
Protein Source:	Cell culture
Concentration:	lot specific
Purity:	Washed harvested organisms are solubilized with detergent to produce a lipoprotein and cytoplasmic cylinder extract enriched with Outer Surface Proteins (OSP)
Buffer:	Presentation State: Lysate State: Liquid Lysate Buffer System: PBS, pH 7.2 containing 1% N-octylglucoside detergent Preservative: None
Preparation:	Liquid Lysate
Applications:	ELISA. Western Blot.
Protein Description:	Agent <i>Borrelia afzelii</i> (PKo)
Note:	Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
Storage:	Upon receipt, store (in aliquots) at -20°C to -80°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)