

## Product datasheet for **AR00029PU-N**

### Cotinine Protein

#### Product data:

Product Type:	Native Proteins
Description:	Cotinine protein, 1 ml
Protein Source:	Synthetic
Purity:	HRP coupled to Cotinine at position 3
Conjugation:	HRP
Buffer:	Presentation State: Liquid State: Liquid HRP Buffer System: 0.05M Tris containing protein stabilizers containing 0.005% Thimerosal and other additives Presentation Label: HRP
Preparation:	Liquid HRP
Protein Description:	Cotinine-3-HRP Conjugate with Horseradish Peroxidase
Note:	Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.
Storage:	Store the protein undiluted at 2-8°C. <b>DO NOT FREEZE!</b>
Stability:	Shelf life: one year from despatch.



[View online »](#)