

Product datasheet for **AR00032PU-N**

Coxsackievirus B1 (Strain Conn) Protein

Product data:

Product Type:	Native Proteins
Description:	Coxsackievirus B1 (Strain Conn) protein, 1 mg
Protein Source:	Vero
Concentration:	lot specific
Purity:	Infected cells are disrupted and clarified to remove cell debris. Virus is concentrated by ultracentrifugation and presented in cell culture medium without FCS.
Buffer:	Presentation State: Purified State: Liquid partially purified protein without preservatives
Preparation:	Liquid partially purified protein without preservatives
Protein Description:	Coxsackievirus B1 Antigen, Strain Conn, ATCC VR-28
Note:	Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
Storage:	Store the protein at -20°C to -80°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)