

## Product datasheet for **AR00033PU-N**

### Coxsackievirus B5 (Strain Faulkener) Protein

#### Product data:

<b>Product Type:</b>	Native Proteins
<b>Description:</b>	Coxsackievirus B5 (Strain Faulkener) protein, 1 mg
<b>Protein Source:</b>	Vero
<b>Concentration:</b>	lot specific
<b>Purity:</b>	Infected cells are disrupted and clarified to remove cell debris. Virus is concentrated by ultracentrifugation and presented in cell culture medium without FCS. The resulting antigen preparation consists of a high concentration of virus and viral components as well as some cellular material.
<b>Buffer:</b>	Presentation State: Purified State: Liquid partially purified protein without preservatives
<b>Preparation:</b>	Liquid partially purified protein without preservatives
<b>Protein Description:</b>	Coxsackievirus B5 Antigen, Strain Faulkener, ATCC VR-185
<b>Note:</b>	Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
<b>Storage:</b>	Store the protein at -20°C to -80°C. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.



[View online »](#)