

Product datasheet for **AR00034PU-N**

Coxsackievirus B6 (Strain Schmitt) Protein

Product data:

| | |
|-----------------------------|---|
| Product Type: | Native Proteins |
| Description: | Coxsackievirus B6 (Strain Schmitt) protein, 1 mg |
| Protein Source: | Vero |
| Concentration: | lot specific |
| Purity: | Infected cells are disrupted and clarified to remove cell debris. Virus is concentrated by ultracentrifugation and presented in cell culture medium without FCS. |
| Buffer: | Presentation State: Purified State: Liquid partially purified protein without preservatives |
| Preparation: | Liquid partially purified protein without preservatives |
| Protein Description: | Coxsackievirus B6 Antigen, Strain Schmitt, ATCC #VR-538`bN |
| Note: | Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product. |
| Storage: | Store the protein at -20°C to -80°C. Avoid repeated freezing and thawing. |
| Stability: | Shelf life: one year from despatch. |



[View online »](#)