

## Product datasheet for **AR00047PU-N**

### Epstein Barr Virus / EBV (p23) Human Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Epstein Barr Virus / EBV p23 recombinant protein, 1 mg
Species:	Human
Expression Host:	E. coli
Concentration:	lot specific
Purity:	>95% pure (10% SDS-PAGE, coomassie staining). GST affinity purification
Buffer:	Presentation State: Purified State: Liquid purified Buffer System: 25mM glycine, pH 9.6 containing 50% glycerol
Preparation:	Liquid purified
Applications:	ELISA. Western blot.
Protein Description:	Recombinant Epstein Barr Virus (EBV) Human Herpes Virus-4 (HHV-4), p23. Contains the HHV-4 p23 regions. Not full length p23 protein. Contains GST fusion partner. Immunoreactive with EBV positive sera.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)