

Product datasheet for **AR00053PU-N**

Feline Immunodeficiency Virus (FIV) Human Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Feline Immunodeficiency Virus (FIV) recombinant protein, 0.1 mg
Species:	Human
Expression Host:	E. coli
Concentration:	lot specific
Purity:	90% pure (SDS-PAGE)
Buffer:	Presentation State: Purified State: Liquid purified Buffer System: 50 mM PBS, pH 7.2
Preparation:	Liquid purified
Applications:	ELISA.
Protein Description:	Recombinant Feline Immunodeficiency Virus (FIV) Core Antigen
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)