

Product datasheet for **AR00063PU-N**

Folate receptor alpha Bovine Protein

Product data:

Product Type:	Native Proteins
Description:	Folate receptor alpha bovine protein, 1 mg
Species:	Bovine
Protein Source:	Milk
Concentration:	lot specific
Purity:	>95% pure (HPLC)
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 120 mM Sodium phosphate, pH 7.4 (+/- 0.1) containing 30 mM Sodium Chloride
Bioactivity:	Specific: 18.6 µg/ml (Folate binding)
Preparation:	Liquid purified protein
Protein Description:	Folate Binding Protein (Bovine)
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
RefSeq:	NP_001193456
Locus ID:	516067
Synonyms:	FBP; FOLR1; FR-alpha
Summary:	folate receptor 3 (gamma)



[View online »](#)