

## Product datasheet for **AR00082PU-N**

### Heterophile Antigen Bovine Protein

#### Product data:

Product Type:	Native Proteins
Description:	Heterophile Antigen bovine protein, 1 mg
Species:	Bovine
Protein Source:	Blood
Concentration:	lot specific
Purity:	Single band by SDS-PAGE (reducing conditions)
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: Deionized water
Preparation:	Liquid purified protein
Applications:	<b>ELISA:</b> Titer > 1/40,000.
Protein Description:	Bovine Heterophile Antigen
Storage:	Upon receipt, store undiluted (in aliquots) at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)