

Product datasheet for **AR00087PU-N**

Mouse IgA (kappa) Mouse Protein

Product data:

| | |
|----------------------|---|
| Product Type: | Native Proteins |
| Description: | Mouse IgA kappa mouse protein, 1 mg |
| Species: | Mouse |
| Protein Source: | Myeloma |
| Concentration: | lot specific |
| Purity: | Salt fractionation, ion exchange and immunoaffinity chromatography |
| Buffer: | Presentation State: Aff - Purified State: Liquid purified protein Buffer System: TBS, pH 7.2 containing 0,09% sodium azide |
| Preparation: | Liquid purified protein |
| Protein Description: | Mouse IgA Kappa. Reacts only with anti IgA and not with anti IgG2b, IgG1, IgG3, IgG2a, IgM or IgD on IEP technique when tested at 1mg/ml. |
| Storage: | Store the protein undiluted at 2-8°C. DO NOT FREEZE! |
| Stability: | Shelf life: one year from despatch. |



[View online »](#)