

Product datasheet for **AR00095PU-N**

Influenza A H3N2 (Strain A/Brisbane/10/06) Protein

Product data:

Product Type:	Native Proteins
Description:	Influenza A H3N2 Strain A/Brisbane/10/06 protein, 1 mg
Protein Source:	Chicken eggs
Purity:	>90% pure by SDS-PAGE (Ultracentrifugation with 10-40% sucrose gradient).
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 0.05M Tris-HCl, pH 8.0 containing 0.1M Sodium Chloride, 5 mM EDTA, 0.09% Sodium Azide and 0.005% Thimerosal
Preparation:	Liquid purified protein
Applications:	Serological studies of Influenza A virus, immunogen for antibody production.
Protein Description:	Influenza A Virus, Strain A/Brisbane/10/06 (H3N2)
Note:	Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)