

## Product datasheet for **AR00096PU-N**

### Influenza A H3N2 (Strain A/Shangdong/9/93) Protein

#### Product data:

<b>Product Type:</b>	Native Proteins
<b>Description:</b>	Influenza A H3N2 Strain A/Shangdong/9/93 protein, 1 mg
<b>Protein Source:</b>	Chicken eggs
<b>Concentration:</b>	lot specific
<b>Purity:</b>	>90% pure (SDS-PAGE). Ultracentrifugation using 10-40% sucrose gradient
<b>Buffer:</b>	Presentation State: Purified State: Liquid purified protein Buffer System: 0.05M Tris-HCl buffer, pH 8.0, 0.1M Sodium chloride, 5 mM EDTA containing 0.09% Sodium azide, 0.005% Thimerosal
<b>Preparation:</b>	Liquid purified protein
<b>Applications:</b>	Serological studies of influenza A virus, immunogen for antibody production. Has been tested with anti-influenza A monoclonal antibodies in ELISA.
<b>Protein Description:</b>	Influenza virus, Type A, Strain A/Shangdong/9/93, (H3N2)
<b>Note:</b>	Caution: This product has been treated in a manner consistent with methods of inactivation. Generally accepted good laboratory practices appropriate to microbiological/viral safe handling practices and techniques are required when handling this product.
<b>Storage:</b>	Store the protein at -20°C. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.



[View online »](#)