

Product datasheet for **AR00097PU-N**

Influenza A H1N1 (A/Solomon Islands/03/06) Protein

Product data:

Product Type:	Native Proteins
Description:	Influenza A H1N1 A/Solomon Islands/03/06 protein, 1 mg
Protein Source:	Chicken eggs
Purity:	>90% pure (SDS-PAGE). Ultracentrifugation with 10-40% sucrose gradient
Buffer:	Presentation State: Purified State: Liquid purified protein Buffer System: 0.05M Tris-HCl, pH 8.0 containing 0.1M Sodium chloride, 5 mM EDTA, 0.09% Sodium azide and 0.005% Thimerosal
Preparation:	Liquid purified protein
Protein Description:	Influenza A Strain A/Solomon Islands/03/06 (H1N1)
Storage:	Store the protein at -20°C. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)