

## Product datasheet for **AR00099PU-N**

### Influenza A H3N2 (Strain A/Wisconsin/67/05) Protein

#### Product data:

<b>Product Type:</b>	Native Proteins
<b>Description:</b>	Influenza A H3N2 Strain A/Wisconsin/67/05 protein, 1 mg
<b>Protein Source:</b>	Chicken eggs
<b>Purity:</b>	>90% pure (SDS-PAGE). Ultracentrifugation with 10-40% sucrose gradient
<b>Buffer:</b>	Presentation State: Purified State: Liquid purified protein Buffer System: 0.05M Tris-HCl, pH 8.0 containing 0.1M Sodium chloride, 5 mM EDTA, 0.09% Sodium azide and 0.005% Thimerosal
<b>Preparation:</b>	Liquid purified protein
<b>Protein Description:</b>	Influenza A Virus, Strain A/Wisconsin/67/05 (H3N2)
<b>Note:</b>	Caution: No test guarantees a product to be non-infectious. All materials should be handled as if potentially infectious. Generally accepted laboratory practices appropriate for infectious materials should be employed when handling this product.
<b>Storage:</b>	Store the protein at -20°C. Avoid repeated freezing and thawing.
<b>Stability:</b>	Shelf life: one year from despatch.



[View online »](#)